## Clean-Trace™ Test

-  (EN) Surface Protein High Sensitivity Test
-  (FR) Test Surface Protein Haute Sensibilité
-  (DE) Oberflächen Protein High Sensitivity Test
-  (IT) Test Surface Protein ad alta sensibilità
-  (ES) Test de alta Sensibilidad para Proteínas en Superficie
-  (NL) Oppervlakte hoge gevoeligheid eiwittest
-  (SV) Surface Protein High Sensitivity Test
-  (DA) Overflade Protein Høj Sensitivitet Test
-  (NO) Høysensitiv proteintest til overflate
-  (FI) Surface Protein Herkkä pintatesti
-  (PT) Teste de Proteína de Alta Sensibilidade de Superfície
-  (EL) Τεστ Επιφανειών Protein Υψηλής Ευαισθησίας
-  (PL) Test o wysokiej czułości do detekcji białka na powierzchni
-  (HU) Felszíni Fehérje Nagy érzékenységű Teszt
-  (CS) Povrchový protein vysoký citlivostní test
-  (RO) Proteine de suprafață Test de înaltă sensibilitate
-  (RU) Тест для определения белка на поверхности (высококчувствительный)
-  (TR) Yüksek Hassasiyetli Yüzey Protein Testi
-  (JA) 高感度タンパク残留測定
-  (ZH) 高灵敏度表面蛋白质采样棒

# Clean-Trace™ Test

# Product Instructions

## Surface Protein High Sensitivity Test

### Product Description and Intended Use

3M™ Clean-Trace™ Surface Protein High Sensitivity Test detects protein residues. The results of the tests are time and temperature dependent. By adjusting the time and / or temperature it is possible to achieve different limits of detection. Selecting the correct incubation time and temperature can therefore optimize the sensitivity of the test. A 45 minute incubation at 37°C (99°F) results in equivalent sensitivity to the Ninhydrin test specified in the UK HTM2030 guidelines. If your requirements for sensitivity are different please use the Time Temperature matrix below as a guide.

Temperature / Time (Minutes) / Sensitivity (micrograms)				
	10 Minutes	15 Minutes	30 Minutes	45 Minutes
20°C (68°F)	50 µg	30 µg	20 µg	10 µg
37°C (99°F)	20 µg	10 µg	5 µg	3 µg

Based on the Biuret reaction. Under alkaline conditions the copper ions (Cu<sup>2+</sup>) form a complex with peptide bonds of protein and become reduced to Cu<sup>+</sup>. Bicinchoninic acid (BCA) under alkaline conditions is a highly sensitive, stable and specific reagent for Cu<sup>+</sup> forming a purple complex. The test provides an estimation of surface protein residue by means of a color change result.

The color may form on the swab more than the solution with certain types of soil. Record the strongest color change on the swab or in the solution as the end result.

3M Clean-Trace Surface Protein High Sensitivity Test can be used to assess the level of contaminating protein residue on instruments after they have been processed in a washer/ disinfectant. The test should be used after cleaning to assess whether instruments are clean enough for effective sterilization.

### ⚠ Warnings

Consult the Safety Data Sheet for additional information and local regulations for disposal.

### ⚠ Cautions

Protein is a common substance. It is present on skin, hair and on many surfaces. Do not touch the swab or the surface to be tested.






Do not use 3M Clean-Trace Surface Protein High Sensitivity Test past expiry date.

Although the 3M Clean-Trace Surface Protein High Sensitivity Test is very sensitive, a surface cannot be considered sterile on the basis of a 3M Clean-Trace Surface Protein-High Sensitivity Test result.

There is minimal interference with most common sanitizers at normal concentrations. However high concentrations of alkaline cleaners may cause false negative results in some instances. Peroxide based disinfectants can cause a false positive color change reaction with the 3M Clean-Trace Surface Protein High Sensitivity Test.

### User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at [www.3M.com/foodsafety](http://www.3M.com/foodsafety), or contact your local 3M representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any 3M Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

### **Limitation of Warranties / Limited Remedy**

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, 3M DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any 3M Food Safety Product is defective, 3M or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify 3M within sixty days of discovery of any suspected defects in a product and return it to 3M. Please call Customer Service (1-800-328-1671 in the U.S.) or your official 3M Food Safety representative for a Returned Goods Authorization.

### **Limitation of 3M Liability**

3M WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall 3M's liability under any legal theory exceed the purchase price of the product alleged to be defective.

### **Storage and Disposal**

For maximum shelf life, store between 2-25°C (36-77°F). Alternatively, swabs are stable for 14 days when shipped or stored at elevated temperature; for this storage option do not exceed 35°C (95°F). Expiry date and lot number are noted on each pouch. After use, may contain microorganisms that may be a potential biohazard. Follow current industry standards for disposal.

### **Instructions for Use**

Before beginning the test procedure, ensure the dry bath, water bath, or incubator is turned on and set for 37°C (99°F).

1. Place the 3M Clean-Trace Surface Protein High Sensitivity Test at room temperature for at least 10 minutes before use. Grip the sample stick handle and remove the swab from the device. Use the moisturizer provided and apply 4 drops to the swab OR apply 4 drops to the surface to be tested [figure 1].
2. Swab the test instrument. Where practical, swab an area approximately 10 cm by 10 cm and swab the area in one direction and then swab the area in the opposite direction [figure 2]. Test hard to reach places like serrated edges and joints. Apply pressure on the swab and rotate the swab as the sample is collected to ensure repeatable and effective sampling [figure 3].
3. Reinsert the sample stick into the test with the handle inserted to the original position of the unused test. At this point you can either activate and measure immediately, or, if it is more convenient, the 3M Clean-Trace Surface Protein High Sensitivity Test may be labeled and left un-activated for up to four hours before activation and measurement of multiple swabs simultaneously.
4. To process the sample, push down firmly on the top of the 3M Clean-Trace Surface Protein High Sensitivity Test sample stick handle [figure 4]. The handle will slide into the device tube and the top of the handle should be level with the top of the device tube when fully depressed. Grip the top of the device and shake rapidly side-to-side for at least five seconds to mix the sample and reagent [figure 5]. The solution will turn mint green.
5. Place the 3M Clean-Trace Surface Protein High Sensitivity Test into a heating block that has been pre-heated to 37°C (99°F). Heat for 45 minutes at 37°C (99°F). Compare the color of the swab and the solution against the 3M Clean-Trace Surface Protein High Sensitivity Test label. Record the strongest color change in the solution or swab as the end result.

How to best utilize 3M Clean-Trace Surface Protein High Sensitivity Test	
Do	Do Not
Hold the 3M Clean-Trace Surface Protein High Sensitivity Test in the VERTICAL position on activation and shake rapidly from SIDE TO SIDE for at least 5 seconds.	Activate when not in the vertical position, do not shake up and down, do not shake aggressively.
Incubate 3M Clean-Trace Surface Protein High Sensitivity Test vertically.	Incubate 3M Clean-Trace Surface Protein High Sensitivity Test horizontally.
Read 3M Clean-Trace Surface Protein High Sensitivity Test at 45 minutes.	Record any color change after 45 minutes.

**Interpretation**

The color of the test reaction will indicate the level of protein residues on the surface. By comparing the color produced against the 3M Clean-Trace Surface Protein High Sensitivity Test label an estimation of the surface cleanliness can be made.

- Green indicates a Pass result - Clean, No further action required
- Grey or Purple indicates a Fail result, Re-clean & Re-test

For further information, contact 3M Food Safety Technical Services.

**Explanation of Symbols**

[www.3M.com/foodsafety/symbols](http://www.3M.com/foodsafety/symbols)

## 3M Food Safety

### 3M United States

3M Center  
Bldg. 275-5W-05  
St. Paul, MN 55144-1000  
USA  
1-800-328-6553

### 3M Canada

Post Office Box 5757  
London, Ontario N6A 4T1  
Canada  
1-800-563-2921

### 3M Latin America

3M Center  
Bldg. 275-5W-05  
St. Paul, MN 55144-1000  
USA  
1-954-340-8263

### 3M Europe and MEA

3M Deutschland GmbH  
Carl-Schurz-Strasse 1  
D41453 Neuss/Germany  
+49-2131-14-3000

### 3M United Kingdom PLC

Morley Street, Loughborough  
Leicestershire  
LE11 1EP  
United Kingdom  
+(44) 1509 611 611

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Euro Plaza  
Gebaude J, A-1120 Wien  
Kranichberggasse 4  
Austria  
+(43) 1 86 686-0

### 3M Asia Pacific

No 1, Yishun Avenue 7  
Singapore, 768923  
65-64508869

### 3M Japan

3M Health Care Limited  
6-7-29, Kita-Shinagawa  
Shinagawa-ku, Tokyo  
141-8684 Japan  
81-570-011-321

### 3M Australia

Bldg A, 1 Rivett Road  
North Ryde, NSW 2113  
Australia  
61 1300 363 878

### 3M Health Care

2510 Conway Ave  
St. Paul, MN 55144 USA  
[www.3M.com/foodsafety](http://www.3M.com/foodsafety)

© 2018, 3M. All rights reserved.  
3M and Clean-Trace are trademarks of 3M. Used under license in Canada.  
34-8721-8853-6

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### 3M Österreich GmbH